Solitaire® 2

DE	Gebruiksaanweisung	SE	Bruksanvisning
GB	Instructions for use	DK	Brugervejledning
FR	Mode d'emploi	NO	Bruksanvisning
ES	Instrucciones de uso	FI	Käyttöohjeet
IT	Istruzioni per l'uso	GR	Οδηγίες Χρήσης
PT	Instruções de uso	CZ	Návod k použití
NL	Gebruiksaanwijzing	HU	Használati utasítás

LV	Lietošanas instrukcija	Manufactured by: Heraeus Kulzer GmbH Grüner Weg 11 63450 Hanau (Germany) info@heraeus.com www.heraeus-dental.com
LT	Naudojimo instrukcija	
PL	Instrukcja obsługi	
TR	Kullanma talimatı	
RU	Инструкция по применению	

Heraeus

Distributed in USA/Canada exclusively by:
Heraeus Kulzer, LLC
300 Haraeus Way
South Bend, IN 46614-2517
1-800-431-1785
Caution: Federal law restricts the sale by or on the order of a dentist.

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Solitaire® 2

Instructions for use	GB
Light curing radiopaque Polyglas® hybrid composite for posterior restorations	
Solitaire® 2 is a lightcuring, radiopaque Polyglas® hybrid composite. The final polishing may easily be realized, it is abrasionproof and in connection with adequate bonding systems it is meant for adhesive filling therapy.	
Shades: A1, A2, A3	
Composition: Solitaire® 2 is produced on the basis of multiply cross-linking urethane (meth)acrylate monomers and contains the following inorganic fillers at 58%(v/v) with a particle size distribution of 0.02–23 µm: • Ba Al-F silicate glass (Ø 0.7 µm; max. <2 µm) • Porous silicon dioxide (Ø 8 µm; max. <23 µm)	
Applications: • Fillings of the classes I and II according to Black • Fillings of class V according to Black • Splinting of teeth which are traumatically loose • Restorations of primary teeth • Reconstruction on stumps	
Not intended for indications or uses that are not cited in these instructions for use.	

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Side effects: Ce produit ou one of its components may in particular cases cause hypersensitive reactions. If suspected, information on the ingredients can be obtained from the manufacturer	
Kind of application: In order to prepare the treatment it is recommended to clean the hard dental substance by means of a polishing paste which is free of fluoride. Prior to starting the treatment the shades should be selected from the attached color scale as long as the teeth are still humid. Standard plugging devices are suitable for application. 1. Preparation of the cavity according to the general regulations of the composite filling therapy. Then dry and clean the cavity. Attention: Take care of a dry field of operation. 2. In case of many-faced cavities use a metal matrix – also for the proximal area – and wedge it carefully in order to avoid cervical excess and to achieve optimal contacts for the points of contact. 3. Prior to applying the filling material to the cavity it is necessary to apply an adequate adhesive system (e.g. iBOND® Self Etch, iBOND Total Etch by Heraeus Kulzer). Avoid contamination of saliva. 4. Apply Solitaire® 2 to the cavity in thin layers of max. 2 mm, condense with a plugging device (e.g. ball plunger) and adapt carefully to the cavity walls. In general, Solitaire® 2 may be applied in horizontal layers. Merely in case of large class II cavities, especially for amalgam replacement cavities with remaining cavity walls of 2 mm or less it is recommended to use a transverse incremental hardening technique. 5. Polymerize Solitaire® 2 in thin layers of max. 2 mm for a period of 40 s using a Translux® light device. Approach recess side of light as close as possible to the filling. During polymerization a dispersive layer is built up at the surface which may not be touched or removed. It serves as connection to the following polyglass layer. After the matrix has been removed it is recommended to separately polymerize the cervical areas of the fillings for 40 s because these areas may lay in the shadow of the matrix. 6. After polymerization Solitaire® 2 may immediately be processed. Finishing diamonds and flexible discs are suitable for processing.	

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Further recommendations and precautions: • When building up large restorations or those, which require lengthy time to complete, the overhead operating light should temporarily be turned off or aimed away from the working area in order to prevent premature curing of the composite. • The light curing times mentioned above are the minimum times recommended and should not be reduced in order to achieve optimal restoration quality and physical properties of the material. • The PLT's (pre-loaded tips) are intended for single use only and cannot be sterilized or disinfected without damaging and contaminating the restorative material because the solutions contain phenolic constituents. • PLT's are for single patient use only. • PLT's should

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